Polylaser trion
Polylaser brush
Polylaser derma

Mode d'emploi

Table des matières

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Radiated disturbance measurements	Conformity	Electromagnetic ambient conditions - guideline
HF-radiation to CICPR 11	Group 1	The laser device uses HF-energy only for its internal function. Therefore HF radiation is very low. The influence on neighbouring electric devices is unlikely...
HF-radiation to CICPR 11	Class B	The laser device may be used in all kinds of buildings, including residential buildings and those directly connected to mains for the public supply which are also servicing buildings used for residential purposes.
Radiation of harmonic oscillation according to EC 61000-3-2	Conforms	
Radiation of voltage variations/ Flicker effect according to IEC 61000-3-3	Not applicable	

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Immunity tests	IEC 60601-Test level	Conformity level	Electromagnetic ambient conditions - guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	conforms	The floor should be made from wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material the relative humidity should be minimum 30%
Fast transient electric disturbance variable / Bursts according to IEC 61000-4-4	± 2kV for network leads ± 1 kV for input and output leads	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions.
Impulse voltage (Surges) according to IEC 61000-4-5	± 1 kV opposed mode voltage ± 2 kV common mode voltage	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions.
Voltage dips, short interruptions and voltage variations according to IEC 61000-4-11	< 5% U_T (> 95% dip of U_T) for ½ period 40% U_T (60% dip of U_T) for 5 periods 70% U_T (30% dip of U_T) for 25 periods < 5% U_T (> 95% dip of U_T) for 5 s	conforms	The quality of the supply voltage should be typical for business or hospital ambient conditions. If the user wants to continue operation in case of an interrupted energy supply it is recommended to feed the laser device from a system without interruption or a battery.
Magnetic field at supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	conforms	Magnetic fields at mains frequency should have the typical values for business and hospital ambient conditions.

Note: U_T .alternating voltage mains before the application of the test level.

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Immunity tests	IEC 60601-test level	Conformity level	Electromagnetic ambient conditions - guidelines
HF-disturbance variables (transmission) according to IEC 61000-4-3	3 Veff 150 kHz up to 80 MHz	conforms	Portable and mobile radio sets must not be operated in a distance to the laser device and the leads less than the recommended protective distance what was calculated according to the equation for the transmission frequency. Recommended protective distance: $d = 1,2 \sqrt{P}$
HF-disturbance variables (radiation) according to IEC 61000-4-3	3 V/m 80 MHz up to 2,5 GHz	conforms	$d = 1,2 \sqrt{P}$ for 80 MHz up to 800 MHz $d = 2,3 \sqrt{P}$ for 800 MHz up to 2,5 GHz with P as power rating of the transmitter in Watt (W) according to data from the manufacturer of the transmitter and as recommended protective distance in meter (m). According to a local ^a investigation the field intensity for all frequencies of immobile radio transmitters should be less than the conformity level ^b . Interference is possible in the surroundings of devices carrying the following signs. 

Note 1 At 80 MHz and 800 MHz the higher frequency range is applied.

Note 2 These guidelines are not applicable in all cases. The propagation of electromagnetic quantities is influenced by the absorption and reflection of buildings, items and human beings.

^a It is not possible to predict theoretically the field strength of immobile transmitters like base stations of mobile phones and mobile land radio sets, amateur radio stations, AM and FM television and radio broadcasting transmitters. In order to determine electromagnetic ambient conditions with regard to immobile transmitters a study of the location should be taken into consideration. If the field strength measured in the place where the laser device is operated exceeds the above mentioned conformity level the laser device should be watched to proof its function in accordance to the requirements. Unusual performance data may require additional measures, for instance change of orientation or different location of the laser device.

^b Above the frequency range from 150 kHz to 80 MHz the field strength should be less than 3 V/m.

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Recommended protective distance between portable and mobile HF-telecommunication devices and the laser device			
The laser device is designed for the operation in electromagnetic ambient conditions where HF-disturbance variables are controlled. Client or user of the laser device can help to avoid electromagnetic disturbances by keeping the minimum distance between portable and mobile HF-telecommunication devices (transmitters) and the laser device as indicated below, dependent on the power output of the communication device.			
Power rating of the transmitter W	Protective distance dependent on the transmitting frequency m		
	150 kHz up to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz up to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz up to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
<p>If the maximum power rating for certain transmitters is not given in the above table the recommended protective distance d in metres (m) can be determined by using the equation indicated in the respective column, where P is the maximum power rating of the transmitter in Watt (W) according to the manufacturer of the transmitter.</p> <p>Note 1 At 80 MHz and 800 MHz the higher frequency range is applied</p> <p>Note 2 These guidelines are not applicable in all cases. The propagation of electromagnetic quantities is influenced by the absorption and reflection of buildings, items and human beings.</p>			

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REIMERS & JANSSEN GmbH